

Determination of Cell Viability by 7-AAD Exclusion Using Flow Cytometry

SOP-QC-0142 · Version 1.0 (DRAFT) · Quality Control – Analytical · GMP / 21 CFR 211, data integrity per ALCOA+

1. Purpose

To describe the procedure for quantifying cell viability in cell-therapy drug product using 7-aminoactinomycin D (7-AAD) exclusion measured by flow cytometry.

2. Scope

Applies to in-process and release viability testing of cell-therapy products by the QC Analytical group. Excludes viability by trypan blue and automated cell counters, which are covered under separate procedures.

3. Responsibilities

- QC Analyst — executes the procedure, records results contemporaneously, and reports out-of-specification results.
- QC Reviewer — performs the second-person review of data and calculations per data-integrity requirements.
- QC Supervisor — approves results and dispositions out-of-specification investigations.

4. Materials & Systems

- Calibrated flow cytometer qualified for the viability application
- 7-AAD viability staining solution (record lot and expiry)
- Stain buffer (PBS + 2% FBS or equivalent), 2–8 °C
- Single-use sterile tubes, pipettes, and tips
- Validated acquisition/analysis software with audit-trail enabled

5. Procedure

- 1 Confirm the flow cytometer has passed daily QC and is within its qualification interval; record instrument ID and QC status.
- 2 Record the 7-AAD reagent lot number and expiry; confirm it is within date.
- 3 Prepare a single-cell suspension and adjust to the target concentration specified in the batch record.
- 4 Transfer the specified cell volume to a labeled tube; prepare an unstained control in parallel.
- 5 Add the validated volume of 7-AAD staining solution to the test tube; do not add stain to the unstained control.
- 6 Incubate protected from light for the validated time at the specified temperature.
- 7 Acquire the unstained control first to set the negative gate, then acquire each test sample, collecting the specified minimum event count.
- 8 Gate the viable (7-AAD-negative) population using the strategy defined in the method; export the % viable result.
- 9 Record the result contemporaneously in the batch record and route data for second-person review.

6. Acceptance Criteria

- Minimum acquired events meet or exceed the method-defined threshold.
- Unstained control confirms correct negative-gate placement before test acquisition.
- Viability result is reported to one decimal place with the associated specification.

7. References

- Method: Flow Cytometry Viability Method (7-AAD) — [method number]
- SOP: Operation and Daily QC of the Flow Cytometer — [SOP number]
- Data Integrity Policy — [policy number]

8. Revision History

Version	Date	Summary of Change
1.0 (DRAFT)	Pending	Initial issue.

QC Review Record

#	Category	Review Item	Status
1	Completeness	All required sections present (Purpose, Scope, Responsibilities, Materials, Procedures, Acceptance, References, Revision History).	PASS
2	Completeness	Document ID, version, and effective-date fields populated.	REVIEW
3	Clarity	Each procedure step is a single, concrete action written in the imperative.	PASS
4	Clarity	No undefined acronyms on first use.	PASS
5	Compliance	Second-person review step is explicitly required in the procedure.	PASS
6	Data Integrity (ALCOA+)	Contemporaneous recording of results is specified.	PASS
7	Data Integrity (ALCOA+)	Reagent lot/expiry and instrument ID/QC status are captured (attributable, original).	PASS
8	Human Decision Required	Confirm validated stain volume, incubation time/temperature, and minimum event count before use.	HUMAN DECISION
9	Human Decision Required	Verify the exact specification/acceptance limit for viability for this product.	HUMAN DECISION
10	References	Referenced method/SOP/policy numbers are placeholders and must be completed.	HUMAN DECISION

Status: DRAFT — NOT approved. Human QC items outstanding.

Reviewed by: Kalif Shear Date: _____

Generated as a first draft by an AI-assisted, human-in-the-loop workflow (SOP-Forge). Not valid until reviewed, completed, and approved by a qualified person under the site quality system.